

Protocol: Immunohistochemical staining (IHC)

MATERIALS/ REAGENTS/ BUFFERS

- Tissue sections
- Primary antibody
- Conjugated secondary antibody
- ABC reagents
- Xylene
- 100% ethanol
- 95% ethanol
- 3% hydrogen peroxide (H₂O₂)
- PBS
- TBS
- Antigen retrieval solution: Ex. Sodium Citrate Solution (10mM Sodium Citrate, 0.05% Tween 20, pH 6.0)
- Pepsin
- Blocking buffer: 5-10% normal animal serum

PROTOCOL

I. Tissue fixation

1. Harvest tissue and rinse in Phosphate-buffered saline (PBS) to remove blood.
2. Cut tissue into 3 mm slices before transfer into formalin due to the slow rate of diffusion of neutral buffered formalin (NBF).
3. Place tissue in at least 20 volumes of fixative 10% NBF.
4. Incubate for necessary fixation time (fixation between 18-24 hr seems to be ideal for most applications).
5. Rinse tissue with PBS, store at 4°C in 75% ETOH in H₂O.

II. Tissue processing

1. After fixation, rinse tissue with PBS until fixative is completely removed.
2. Dehydrate tissue using ethanol in the following sequence:
 - 50% Ethanol: 50 min
 - 75% Ethanol: 50 min
 - 90% Ethanol: 50 min
 - 95% Ethanol: 50 min
 - 95% Ethanol: 90 min
 - 100% Ethanol: twice, each for 50 min
 - 100% Ethanol: 90 min
3. Exchange ethanol with xylene in the following sequence
 - Ethanol : Xylene (1:1): 50 min
 - Xylene : 50 min
 - Xylene : 50 min

II. Tissue processing

3. Exchange ethanol with xylene in the following sequence
Ethanol : Xylene (1:1): 50 min
Xylene : 50 min
Xylene : 50 min
4. Exchange xylene with paraffin at 58°C.
Xylene : Paraffin (1:1): 60 min
Paraffin: 120 min
Paraffin: 120 min

III. Tissue embedding

1. Put small amount of molten paraffin in mold, dispensing from paraffin reservoir
2. Transfer tissue into mold by warm forceps, placing cut side down.
3. Transfer mold to cold plate, paraffin will solidify in a thin layer which holds the tissue in position.
4. Add hot paraffin to the mold from the paraffin dispenser. Be sure there is enough paraffin to cover the face of the plastic cassette.
5. Paraffin should solidify in 30 minutes.

IV. Tissue Sectioning

1. Placed tissue blocks face down on an ice block or cold plate for 10-30 minutes.
2. Place a fresh blade on the microtome. Blades may be used to section up to 10 blocks, but replace if sectioning becomes problematic
3. Insert the block into the microtome chuck so the wax block faces the blade and is aligned in the vertical plane.
4. Set the dial to cut 10µm sections to order to plane the block; once it is cutting smoothly, set to 4µm sections .
5. Pick sections and float them on the surface of the cold water. Float the sections onto the surface of clean glass slides.
6. Place the slides with paraffin sections on the surface of the 56°C water bath to bond the tissue to the glass.

V. Immunohistochemical staining (IHC)

Dewaxing & Rehydrate

1. Place the slides with 4µm paraffin sections on the warming block in a 60°C oven for 30 minutes.
2. Perform the following washes
Xylene: twice, each for 5 min
100% ethanol: twice, each for 5 min
95% ethanol: 5 min
75 % ethanol: 5 min
50 % ethanol: 5 min
Rinse slides with deionized water

Antigen retrieval (Heat-Induced Epitope Retrieval)

1. Fill individual plastic container(s) with appropriate retrieval solution and add enough deionized water to pressure cooker chamber or vegetable steamer.
2. Bring slides to a boil in antigen retrieval solution for 15 min in the pressure cooker or 30 min in the vegetable steamer.
3. Transfer slides to room temperature Tris-Buffered Saline (TBS).
4. Tissues are now ready for blocking and primary antibody incubation.



www.genetex.com

International

Tel 886.3.6208988
Fax 886.3.6209098
Address 6F-2, No.89, Dongmei Rd., Hsinchu 300, Taiwan (R.O.C.)
E-mail info@genetex.com

USA

Toll-free 1.877.GeneTex(1.877.436.3839)
Fax 1.949.309.2888
Address 2456 Alton Parkway, Irvine, CA 92606 USA
E-mail info@genetex.com

IHC staining

1. Wash sections in TBS twice, each for 3 min.
2. If using an HRP conjugate for detection, blocking of endogenous peroxidase can be performed here by incubation of sections in 3% H₂O₂ for 10 min.
3. Wash section in TBS twice, each for 3 min.
4. Block each section with blocking buffer for 30 min at RT.
5. Remove blocking buffer and add primary antibody at proper dilution and incubate for overnight at 4°C.
6. Remove antibody solution and wash sections in TBS three times, each for 5 min.
7. Add secondary antibody, and incubate for 30 min at RT.
8. Prepare avidin-biotin-enzyme complex (ABC) reagent according to the manufacturer's instructions.
9. Remove secondary antibody solution and wash sections in TBS three times, each for 5 min.
10. Add ABC reagent to each section and incubate for 30 min at RT.
11. Remove ABC reagent and wash sections in TBS three times, each for 5 min.
12. Add DAB to each section and monitor staining closely.
13. As soon as the sections develop, immerse slides in ddH₂O.
14. Counter-stain sections in hematoxylin.
15. Wash sections in ddH₂O twice, each for 3 min.
16. Dehydrate sections
 - 75 % ethanol: 3 min
 - 95% ethanol: 3 min
 - 100% ethanol: 5 min
 - Xylene: 5 min
17. Mount sections with coverslip.



International

Tel 886.3.6208988
Fax 886.3.6209098
Address 6F-2, No.89, Dongmei Rd., Hsinchu 300, Taiwan (R.O.C.)
E-mail info@genetex.com

USA

Toll-free 1.877.GeneTex(1.877.436.3839)
Fax 1.949.309.2888
Address 2456 Alton Parkway, Irvine, CA 92606 USA
E-mail info@genetex.com

ABOUT




GeneTex

GeneTex started as a small scientist-founded company in San Antonio, Texas, and has grown to become one of the top antibody manufacturers.

Our corporate mission, to provide our customers with quality reagents and to accelerate life sciences research, reflects the philosophy and approach we employ when we manufacture our products. Through extensive research, development, and quality-assurance testing, we have produced and validated a comprehensive collection of antibodies and research reagents.

Due to GeneTex's experience in antibody production, we are uniquely capable of providing technical support and troubleshooting assistance for our customers.

Our Antibodies



GeneTex

www.genetex.com

International

Tel 886.3.6208988

Fax 886.3.6209098

Address 6F-2, No.89, Dongmei Rd., Hsinchu 300, Taiwan (R.O.C.)

E-mail info@genetex.com

USA

Toll-free 1.877.GeneTex(1.877.436.3839)

Fax 1.949.309.2888

Address 2456 Alton Parkway, Irvine, CA 92606 USA


E-mail info@genetex.com

Selection you need

Comprehensive Coverage

Genetex currently offers more than 55,000 products including over 45,000 primary antibodies that cover eleven different research areas. We continue to build upon our product portfolio to ensure that novel targets are available to researchers and that Genetex keeps pace with new scientific discovery.

In addition to our comprehensive catalog of primary antibodies, GeneTex also offers a wide variety of quality reagents, kits and services.



www.genetex.com

International

Tel 886.3.6208988
 Fax 886.3.6209098
 Address 6F-2, No.89, Dongmei Rd., Hsinchu 300, Taiwan (R.O.C.)
 E-mail info@genetex.com

USA

Toll-free 1.877.GeneTex(1.877.436.3839)
 Fax 1.949.309.2888
 Address 2456 Alton Parkway, Irvine, CA 92606 USA
 E-mail info@genetex.com