

Verotoxin / Shiga toxin (SLT-2) subunit B antibody [7582]

Cat. No. GTX02558

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Neutralizing/Inhibition
Reactivity	Bacteria

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:20-1:200
Neutralizing/Inhibition	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified IgG1
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).