

## HIV1 p24 antibody [BDI690]

## Cat. No. GTX41582

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	ELISA, Lateral Flow, Sandwich ELISA
<b>Reactivity</b>	Human Immunodeficiency virus

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX41582, Detection : GTX41583 / GTX40080. The detection antibody requires HRP conjugation. Please contact us for the conjugation service.**

Not tested in other applications.

**Product Note** Specific for HIV-1 p24.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	3.80 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).