

E. coli deoC protein, His tag

Cat. No. GTX57500-pro**Species** E. coli**Package**

100 µg

Applications

Observed MW (kDa) 28 kDa.


Properties

Form Liquid**Buffer** 20mM Tris-HCl, 10% Glycerol, 2mM DTT**Preservative** No preservative**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**
Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MTDLKASSLR ALKLMDLTL NDDDTDEKVI
ALCHQAKTPV GNTAAICIYP RFIPIARKTL KEQGTPEIRI ATVTNFPHGN DDIDIALAET RAAIAYGADE VDVVFPYRAL MAGNEQVGF
LVKACKEACA AANVLLKVII ETGELKDEAL IRKASEISIK AGADFIKTST GKAVAVNATPE SARIMMEVIR DMGVEKTVGF KPAGGVRTAE
DAQKYLAID ELFGADWADA RHYRFGASSL LASLLKALGH GDGKSASSY**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX57500-pro Image**

3 μ g deoC protein (GTX57500-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).