

Zebrafish


Insights into Biology

Quality Antibodies · Quality Results


Neuroscience




Pax2a
Pax2a antibody (GTx128127)




Elavl3+4
Elavl3+4 antibody (GTx128365)




Sox2
Sox2 antibody (GTx124477)




Sox10
Sox10 antibody (GTx128374)



ALDH1A2
ALDH1A2 antibody [C1C3] (GTx115148)




Dvl2
Dishevelled 2 antibody [N2C2] Internal (GTx103878)



GFAP
Gfap antibody (GTx128645)




Hematology & Cardiovascular Research




Hbae1
Hbae1 antibody (GTx128636)




Mpx
Mpx antibody (GTx128379)



Lyz
Lyz antibody (GTx132488)




Mfap4
Mfap4 antibody (GTx132692)




Myl7
Myl7 antibody (GTx128346)


Muscle Research




Acta2
Acta2 antibody (GTx124505)




Myl2
Myl2 antibody (GTx128816)



Tpm4
Tropomyosin 4 antibody (GTx114851)



Dnm2
Dnm2 antibody (GTx127330)




Tnnc1a+b
Tnnc1a+b antibody (GTx127369)


Cell Biology




PcnA
Pcna antibody (GTx124496)




Tp63
Tp63 antibody (GTx124660)




Cdh1
Cdh1 antibody (GTx125890)




Cdh2
Cdh2 antibody (GTx125885)




Ctnnb1
beta Catenin antibody [N1N2-2],N-term (GTx101435)




Histone H2A.XS139ph
Histone H2A.XS139ph(phospho Ser139) antibody [HL1689] (GTx637288)



p53
p53 antibody (GTx128135)



mTOR
Mtor antibody (GTx124771)



Hif1ab
Hif1ab antibody (GTx131826)



Zebrafish



International

Tel 886.3.6208988
Fax 886.3.6208989
Address 6F-2, No.89, Dongmei Rd., Hsinchu 300, Taiwan (R.O.C.)
E-mail international@genetex.com

USA

Toll-free 1.877.GeneTex(1.877.436.3839)
Fax 1.949.309.2888
Address 2456 Alton Parkway, Irvine, CA 92606 USA
E-mail info@genetex.com



GeneTex
www.genetex.com

Hot Products

Cat. No.	Product Name	Clonality	Reactivity	Applications
GT1002199	1α-3-β theta antibody	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P
GT1002637	ACAT1 antibody [N19S]	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P
GT1014233	Acinase 2 antibody [N19S]	Rb,pAb	Hu,Ms,Zfsh,Chk	WB, ICC/F, IHC-P, IP
GT1024484	Acta1b antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Wm, IHC
GT1024455	Acta2 antibody	Rb,pAb	Hu,Zfsh	WB, IHC-P, IHC-Fr, IHC-Wm
GT1024462	Actl1a antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Wm
GT1024770	Adenovirus Hexon 1 antibody [C2C3, C-term]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1027293	Akt3 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1004139	ALAS-E antibody [N3C3]	Rb,pAb	Hu,Ms,Zfsh	WB
GT1028399	Alcama antibody	Rb,pAb	Zfsh	IHC-Wm
GT1023973	ALDH1A1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Cat	WB, ICC/F, IHC-P, IHC-Wm
GT1024349	ALDH2 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC
GT1001029	ALDH2 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm, ELISA, Sandwich, ELISA
GT1014914	Aldolase B antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1011095	alpha 1 Catenin antibody [N3C2], Internal	Rb,pAb	Hu,Ms,Rat,Zfsh,Cat,Dog	WB, ICC/F, IHC-P, IHC-Fr
GT1010076	Alpha 2 Antiplasmin antibody [N3C3]	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1032040	alpha Actinin 1 antibody [N2N5]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1030219	alpha Actinin 2 antibody [N1N3]	Rb,pAb	Hu,Ms,Zfsh,Pig	WB, ICC/F
GT1010160	alpha Actinin antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1000034	alpha Smooth Muscle Actin antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1062802	alpha Tubulin antibody [GT114]	Ms,mAb	Hu,Ms,Rat,Zfsh,Dm,Cruclan carp	WB, ICC/F, IHC-P, IHC-FACS
GT1024303	alpha Tubulin antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Wm, IHC
GT1001984	ARPC2 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1015049	Arp3 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC
GT1014275	ATP synthase gamma antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC
GT1029736	Atp6v1d2 antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1013020	Atp6v1d2 antibody	Rb,pAb	Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1010778	ATP6V1H antibody [N3C3]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1023138	Atp9a1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1023277	Aurkb antibody	Rb,pAb	Zfsh	WB, ICC/F
GT1024588	Bact1n1 antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Fr, IHC-Wm
GT1000064	Bcl-2 antibody [N1N2], N-term	Rb,pAb	Hu,Ms,Rat,Zfsh,Cat,Dog	WB, IHC-P, IHC-Fr, IHC-Wm
GT1063765	beta Actin antibody [HL1926]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh,Dm,Cat,Dog	WB, ICC/F, IHC-P
GT1062930	beta Actin antibody [GT5152]	Ms,mAb	Hu,Ms,Rat,Zfsh,West,Fs,Goat,Dm,Cat,Dog,Mk,Xenopus	WB, ICC/F, IHC-P, IP
GT1001433	beta Catenin antibody [N1N2-2], N-term	Rb,pAb	Hu,Ms,Rat,Zfsh,Rb,Cat,Dog	WB, ICC/F, IHC-P, IHC-Wm, FACS, IP, CHIP assay, IHC
GT1001279	beta Tubulin antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Hm	WB, ICC/F, IHC-P, IHC-Wm
GT1024348	Bmp4 antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC
GT1028200	BmpR1b antibody	Rb,pAb	Zfsh	IHC-Wm
GT1014300	Calmodulin antibody	Rb,pAb	Ms,Rat,Zfsh	WB, IHC-P, IHC
GT1023050	Calpain 1 antibody	Rec,Rb,mAb	Hu,Zfsh	WB, ICC/F, IHC-Fr, FACS
GT1013401	CaMKII delta antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Rb	WB, ICC/F, IHC-P, IP
GT1005594	Caspin Inhibitor 1 alpha 11 antibody	Rb,pAb	Hu,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1035147	Cnna2 antibody	Rb,pAb	Zfsh	WB
GT1033053	Cncl1 antibody	Rb,pAb	Zfsh	WB
GT1025890	Cdh1 antibody	Rb,pAb	Zfsh	IHC-Fr, IHC-Wm, IHC
GT1025885	Cdh2 antibody	Rb,pAb	Zfsh	WB, IHC-Fr, IHC-Wm
GT1028752	Cdh3 antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1000212	CDMP antibody [C2C3], C-term	Rb,pAb	Hu,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1062814	Citrate synthase antibody [G17161]	Ms,mAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1010624	Citrate synthase antibody [N2C3]	Rb,pAb	Hu,Ms,Rat,Zfsh,Chk,Mk,Frg	WB, ICC/F, IHC-P, IHC-Wm, IP
GT1033063	Col1a1a antibody	Rb,pAb	Zfsh	WB
GT1002996	COL1A2 antibody [C2C3], C-term	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC-Wm
GT1027988	Col2a1a antibody	Rb,pAb	Zfsh	IHC-Fr, IHC-Wm
GT1021027	Col3a1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1025882	Crap2a antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1001760	Creatine kinase [brain] antibody [N3C3]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1008969	Creatine kinase [muscle] antibody [N2C3]	Rb,pAb	Hu,Ms,Zfsh,Dog	WB, ICC/F, IHC-P, IHC-Wm, IHC
GT1003627	Creatine kinase MT antibody [N2C3]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1028677	Ctfr antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Fr
GT1014349	CSRP1 antibody [N1C3]	Rb,pAb	Hu,Zfsh	WB, IHC-P, IHC
GT1011246	CSRP1 antibody [N1C1]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-Wm
GT1024378	Ctnna antibody	Rb,pAb	Zfsh	WB
GT1028324	Ctstb antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1013876	Cullin-4a antibody [N1C1]	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC-Wm
GT1023244	Cxcr4b antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Wm
GT1024367	Cyfl2 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P
GT1028377	Cyfla antibody	Rb,pAb	Zfsh	WB
GT1028133	Cytp26a1 antibody	Rb,pAb	Zfsh	WB
GT1025985	Cytp26c1 antibody	Rb,pAb	Zfsh	WB
GT1030225	Cytohesin 1 antibody [N1C1]	Rb,pAb	Hu,Zfsh	WB, IHC-Wm
GT1020983	Cytokeratin 13 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1011298	Cytokeratin 13 antibody [N3C3]	Rb,pAb	Hu,Ms,Zfsh,Dog	WB, ICC/F, IHC-P, IHC-Wm
GT1010311	Cytokeratin 8 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1028277	DeltaD antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028289	Dharma antibody	Rb,pAb	Zfsh	IHC-Wm
GT1030878	Deve1led-2 antibody [N2C2], Internal	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC-Fr, IHC
GT1002509	Dishevelled 3 antibody [N1N2], N-term	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Fr, IP, IHC
GT1012555	Disf1 antibody [N3C3]	Rb,pAb	Hu,Rat,Zfsh	WB, ICC/F, IHC-P
GT1028344	Dixa2 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1063696	DNA polymerase beta antibody [HL1512]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1001947	DNAH8 antibody [N3C3]	Rb,pAb	Hu,Zfsh	WB, ICC/F, IHC-P
GT1021330	Dnm2 antibody	Rb,pAb	Zfsh	WB, ICC/F, IHC-Wm
GT1033152	Dopamine Transporter antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1014610	DSPI antibody [N1C3]	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC
GT1000443	E-Cadherin antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Dog	WB, ICC/F, IHC-P, IHC-Wm, IP
GT1020236	EF1A2 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P
GT1024538	Efta antibody	Rb,pAb	Zfsh	WB, IHC-P
GT1029212	EGR2 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC-Fr
GT1028131	Egr2b antibody	Rb,pAb	Zfsh	WB
GT1024468	EF2a antibody	Rb,pAb	Zfsh	WB, IHC-P
GT10637408	ERF antibody [HL1758]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh	WB, IHC-P
GT1028365	Elav3-4 antibody	Rb,pAb	Hu,Zfsh	ICC/F, IHC-Fr, IHC-Wm
GT1028301	Eng2a antibody	Rb,pAb	Zfsh	IHC-Wm
GT1001003	ENOC1 antibody [N3C3]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm, Dot
GT1013429	ENOC3 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1023664	ENP gamma antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1014392	ETFA antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC
GT1028647	Etv4 antibody	Rb,pAb	Zfsh	WB
GT1028318	Eve1 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1011709	Ezrin antibody [N2C2], Internal	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1028236	Fgf8 antibody	Rb,pAb	Zfsh	IHC-Fr, IHC-Wm
GT1024074	Fhl2 antibody	Rb,pAb	Ms,Zfsh	WB, IHC-P, IHC-Wm
GT1013684	Fibrinolytic antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1028304	Fih antibody	Rb,pAb	Zfsh	IHC-Wm
GT1025996	Fmr1 antibody	Rb,pAb	Zfsh	WB

Hot Products

Cat. No.	Product Name	Clonality	Reactivity	Applications
GT1028273	Foxa2 antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1018867	FXR1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1001189	FYN antibody [N1N3]	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-Wm
GT1001006	GABARAP2 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm
GT1025728	Gabra1 antibody	Rb,pAb	Zfsh	WB, IHC-P
GT1024355	Gabra2 antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Wm
GT1063647	GAD65 + GAD67 antibody [HL1181]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh,Cat,Dog	WB, ICC/F, IHC-P
GT1063693	GAD65 antibody [HL1187]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh,Cat,Dog	WB, ICC/F, IHC-P
GT1001881	GAD7 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-Fr, IHC-Wm, IHC-Free/Freezing
GT1012444	gamma Adaptin antibody [N1N3]	Rb,pAb	Zfsh	WB, IHC-Wm
GT1013286	gamma Tubulin antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm, IP, IHC (Formalin), IHC
GT1063766	GAPDH antibody [HL2062]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh,Cat,Dog	WB, ICC/F, IHC-P, IHC-Wm
GT10627408	GAPDH antibody [GT7299]	Ms,mAb	Hu,Ms,Rat,Zfsh,West,Fs,Goat,Dm,Cat,Dog,Mk,Xenopus	WB, ICC/F, IHC-P, EMSA
GT1024502	Gapdh antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Dm,Cat,Dog,E. coli	WB, ICC/F, IHC-P
GT1028333	Gata1a antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1028174	Gliap antibody	Rb,pAb	Zfsh,Fish	WB, ICC/F, IHC-Fr, IHC-Wm
GT1028280	Gliz2 antibody	Rb,pAb	Zfsh,Xenopus	WB, IHC-Wm
GT1024344	Gli3 antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC
GT1009121	Glutamine synthetase antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm, IP, IHC
GT1014442	GNB1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm, IP, IHC
GT1013203	GPI antibody	Rb,pAb	Hu,Zfsh	WB, ICC/F, IHC-P, IHC-Fr
GT10635477	GPX2 antibody [HL125]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh,Dog	WB, IHC-P
GT1015211	GPX1 antibody [N3C3]	Rb,pAb	Zfsh	WB, ICC/F, IHC-Wm
GT1000580	Grp78 antibody [N2C1], Internal	Rb,pAb	Hu,Ms,Rat,Zfsh,Mosq	WB, IHC-P, IHC
GT1033751	Grp antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028636	Hbae1 antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1000513	HDCA1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm, IP, CHIP assay
GT1031303	HDC13 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm, IP, CHIP assay
GT1023238	He6 antibody	Rb,pAb	Zfsh	WB
GT1005348	Hes6 antibody	Rb,pAb	Zfsh	WB, ICC/F
GT1031826	Hf1ab antibody	Rb,pAb	Zfsh	WB
GT1003707	Hf2 alpha antibody [N1N3]	Rb,pAb	Hu,Zfsh	WB, CHIP assay
GT1027343	Histone H2A.X antibody	Rb,pAb	Zfsh,Neohabranchius furzeri	WB, IHC-Wm
GT10637288	Histone H2A.X39ph (phospho Ser139) antibody [HL168]	Rec,Rb,mAb	Zfsh	IHC-Wm
GT1027342	Histone H2A.X5339ph (phospho Ser139) antibody	Rb,pAb	Hu,Zfsh	WB, ICC/F, IHC-Fr, IHC-Wm, IP, CHIP assay
GT1008273	Histone H2B antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Cat,Dog	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm, IHC
GT1015955	Histone H2B antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC
GT1021184	Histone H3K27me3 (Tri-methyl Lys27) antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm, Dot, EM
GT1028954	Histone H3K4me3 (Tri-methyl Lys4) antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Dm	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm, Dot, CHIP assay, IHC
GT10637402	Histone H3S10ph (phospho Ser10) antibody [HL1752]	Rec,Rb,mAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1028116	Histone H3S10ph (phospho Ser10) antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Bov	WB, ICC/F, IHC-P, IHC-Wm, IP, CHIP assay
GT1025660	Histone H4 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1024375	Hic2 antibody	Rb,pAb	Zfsh	WB
GT1001277	HMG81 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Pig	WB, ICC/F, IHC-P, IHC-Wm, IHC-Free/Freezing
GT1001852	hRNPA1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm, IP, CHIP assay
GT1028322	Hoxb1 antibody	Rb,pAb	Zfsh	IHC-Wm
GT10148	HRP antibody	Rb,pAb	Hu,Rat,Zfsh	WB, ICC/F, IHC-P
GT1007715	HS1BP1 antibody [C2C3], C-term	Rb,pAb	Hu,Zfsh	WB, ICC/F, IHC-P
GT1001144	HSC70 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1010089	HSP60 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh,Dm,Hm	WB, ICC/F, IHC-P, ELISA
GT1012963	Hsp70 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1028403	Hsp70 antibody	Rb,pAb	Zfsh	WB
GT1001448	Hsp90 beta antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm, IP, ELISA-AG
GT1063219	Hsp90 antibody [GT10312]	Rec,Rb,mAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Free/Freezing
GT1015640	Irf6 antibody [N2C1], Internal	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1001232	Irf1 antibody	Rb,pAb	Zfsh	WB
GT1019175	Irf212 antibody [C1C2], Internal	Rb,pAb	Hu,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1001691	ILK antibody [N1C1]	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, IHC-P, IHC-Wm, IP, IHC
GT1008408	Inhibin beta B antibody [C3], C-term	Rb,pAb	Hu,Ms,Zfsh	WB, IHC-P, IHC-Fr
GT1028490	Inr antibody	Rb,pAb	Zfsh	IHC-Wm
GT1001136	Insulin Receptor antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P
GT1028001	Irf1 antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1001067	Irf1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Fr, IHC-Wm
GT1013927	Jab1 antibody [N1C1]	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1031944	Jkb1 antibody	Rb,pAb	Zfsh	WB, IHC-P, IHC-Wm
GT1028420	Jkb2 antibody	Rb,pAb	Zfsh	IHC-P, IHC-Fr, IHC-Wm, IHC
GT1031406	Leptin b antibody	Rb,pAb	Zfsh	WB
GT1017114	LRN7 antibody [N1C3]	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1014532	LOXL1 antibody	Rb,pAb	Hu,Zfsh	WB, IHC-Wm
GT1028244	Lox3a antibody	Rb,pAb	Zfsh	IHC-Wm
GT1032379	Lyc antibody	Rb,pAb	Zfsh	WB, IHC-Wm
GT1028129	Lys antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028514	Mag3a2 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1025115	Mapk3 antibody	Rb,pAb	Zfsh	WB, IHC-P
GT1009740	MCM4 antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IP, CHIP assay
GT1032692	Mfp4 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028382	Mfp antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028363	Mfr antibody	Rb,pAb	Zfsh,Carp	WB, IHC-Wm
GT1028198	Mmp14 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028379	Mpx antibody	Rb,pAb	Zfsh,Allanct salmon	IHC-Fr, IHC-Wm, IHC
GT1013661	MSH6 antibody [C1C2], Internal	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-P, IP, EMSA, IHC
GT1032039	Mt1+Mt2 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1028851	Mtmr2 antibody	Rb,pAb	Zfsh	IHC-Wm
GT1024771	Mtor antibody	Rb,pAb	Hu,Ms,Zfsh	WB, ICC/F, IHC-Wm
GT1014808	Munc18-1 antibody	Rb,pAb	Hu,Ms,Rat,Zfsh	WB, ICC/F, IHC-P, IHC-Wm
GT1028816	Myf2 antibody	Rb,pAb	Zfsh	WB, IHC-Fr, IHC-Wm, IHC
GT1028346	Myf7 antibody	Rb,pAb	Zfsh,Japanese Medaka	WB, IHC-Fr, IHC-Wm
GT1028129	MyoD1 antibody	Rb,pAb	Zfsh	WB
GT1033491	Nang antibody	Rb,pAb	Zfsh	IHC